

# FRESHPOINT GRO-575

# FIVE STAGE REVERSE OSMOSIS SYSTEMS



GRO-575M, (with filter life monitor) shown

## **APPLICATIONS**

Drinking water

## **FEATURES • BENEFITS**

- Reduces more than 90% of contaminants found in water
- Effortless click-in cartridge change-outs
- Color-coded cartridges help maintain proper filter sequence
- 3/8" push fittings for easy installation
- High efficiency membrane technology delivers up to 75 gallons per day

### **MODELS**

- ◆ GRO-575B
- GRO-575M (with TDS monitor)

#### **INSTALLATION TIPS**

- Choose a mounting location suitable to support the full weight of the system when operating.
- Install vertically with the cartridge hanging down.
- System may be installed with cartridge flush to the baseboard if desired as the mechanical assist lever pushes cartridges straight out.


#### **OPERATION TIPS**


- Replace pre/post filters every year and RO membrane every 18 months.
- Place replacement cartridges in proper sequence utilizing color-coded shapes.
- System designed with integrated check valves so cold water line does not need to be shut-off to perform cartridge change-outs

#### **INCLUDES**

- Air gap side faucet
- Installation kit
- Installation and operation manual
- 4.4 gallon storage tank

# FRESHPOINT GRO-575 FIVE STAGE REVERSE OSMOSIS SYSTEMS





#### **SPECIFICATIONS**

- Systems (Model/Part No.) GRO-575B/161112 GRO-575M/161152
- Inlet Size 3/8" QC
- System Dimensions
   17.82" W x 5.46" D x 12.54" H
   (453 mm x 139 mm x 319 mm)
- System Weight 10.70 lbs
- ◆ Contaminant Reduction
  Chlorine Taste and Odor,
  Nominal Particulate, Class I,
  Atrazine, Lead, Lindane,
  Chloroform, Arsenic
  [Pentavalent] <= 50 ppb,
  Barium, Cadmium, Chromium
  [Hexavalent], Chromium
  [Trivalent], Copper, Fluoride,
  Radium 226/228, Selenium,
  Turbidity, TDS, VOC'S,
  Cyst: NSF/ANSI Standards 53
  and 58 certified to reduce cysts
  such as Cryptosporidium and
  Giardia by mechanical means.
- TDS Rejection 96.30%

EPA EST. NO. 082989-CHN-001

- Approximate Filter Life\*
   Pre/Post Filters: 1 year
   RO Membrane: 18 months
- System Production Rate 21.08 gpd (79.77 Lpd)
- Recovery Rating 41.05%
- Certifications
   NSF/ANSI Standards 42/53/58
   CSA Standard B483.1
- ◆ Feed Pressure Range
   40 100 psi (2.75 6.89 bar)
- ◆ Temperature Range 40 - 100°F (4.4 - 37.8°C)
- Replacement Cartridges
   (Model/Part No.)
   F1GC-RC/655117-96
   FDF1-RC/655123-96
   F2B2-RC2/655126-96 (twin pack)
   GR075-RC/655122-96
- \* Filter life dependent on gallon usage and overall incoming water quality.



The Models GRO-575B and GRO-575M are Tested and Certified by NSF International against NSF/ANSI Standards 42, 53, 58 and CSA B483.1 for the reduction of substances specified in the service manual or at www.nsf.org.



13845 Bishops Dr. | Suite 200 | Brookfield, WI 53005 | United States P: 262.238.4400 | Customer Service: 800.279.9404 | tech-support@pentair.com | pentair.com

§For a detailed list of where Pentair trademarks are registered, please visit www.pentair.com/registrations. Pentair trademarks and logos are owned by Pentair plc or its affiliates. Third party registered and unregistered trademarks and logos are the property of their respective owners.